

Determination of Hydroxyl Value

Principle:

Hydroxyl value is the number of milligrams of potassium hydroxide equivalent to the acetic anhydride consumed by one gram of sample. The sample is acetylated with acetic anhydride; the excess acetic anhydride is converted to acetic acid by boiling with water, and then is titrated with aqueous KOH.

Apparatus and Reagents:

1. Acetylating Reagent: Acetic anhydride/Pyridine = 1/3 (V/V). **Prepare freshly.**

CAUTION: See cautions below.

2. 0.7 N Alcoholic KOH
3. 3A Alcohol
4. Phenolphthalein indicator (1.0% in 3A Alcohol)
5. 250 mL Round bottom flask
6. 50 mL Burette
7. Air condenser, 4 ft.

Procedure:

1. Into a 250 mL round bottom flask with ground neck, weigh a sample (see *Sample Size* below) to the nearest 0.1 mg. Use another 250 mL RB flask with ground neck as the blank flask.
2. Pipette 5 mL of acetylating mixture into the sample flask, lubricate the neck with silicone grease and connect both flasks to ground glass tip air condensers.

3. Swirl the flasks thoroughly and put the flasks on a steam bath (previously heated to 100°C). Maintain the bath temperature between 100°C - 120°C for 1 hour. At the end of 1 hour, add 10 mL D.I. water from the top of condenser and continue heating for 10 min.
4. Allow the flasks to cool below 50°C and rinse from the top of condensers with 25 mL 3A alcohol (previously neutralized with 0.1 KOH to phenolphthalein end point). Remove condensers from the flasks and allow the flasks to cool further to room temperature.
5. If the solution is cloudy, or if an insoluble residue is present, add 50 mL of a 50/50 (V/V) mixture of toluene and 3A alcohol.

CAUTION: Toluene is flammable and toxic. Wear safety glasses and gloves.

6. Add 1 mL of phenolphthalein indicator solution and titrate to a pink endpoint with 0.7 N alcoholic KOH.

CAUTION: Phenolphthalein is a carcinogen and irritant. Handle with care. Wear safety glasses and gloves.

Lubrizol Advanced Materials, Inc. / 9911 Brecksville Road, Cleveland, Ohio 44141-3247 / TEL: 800.379.5389 or 216.447.5000

The information contained herein is being furnished for informational purposes only, upon the express condition that the User makes its own assessment of the appropriate use of such information. While the information contained herein is believed to be reliable, no representations, guarantees or warranties of any kind are made as to its accuracy, suitability for a particular application or the results to be obtained herefrom. Lubrizol Advanced Materials, Inc. ("Lubrizol") cannot guarantee how any products associated with this information will perform in

combination with other substances or in the User's process. Due to variations in methods, conditions and equipment used commercially in processing these materials, no warranties or guarantees are made as to the suitability of the information or products for the applications disclosed. Lubrizol shall not be liable and the User assumes all risk and responsibility for any use or handling of any material beyond Lubrizol's direct control. LUBRIZOL MAKES NO WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO,

THE IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. It is the User's sole responsibility to determine if there are any issues relating to patent infringement of any component or combination of components relating to the supplied information. **Nothing contained herein is to be considered as permission, recommendation, nor as an inducement to practice any patented invention without permission of the patent owner.**

For further information, please visit: www.lubrizol.com/personalcare

Calculations:

$$\text{Hydroxyl Value} = \frac{[\text{mL}(\text{blank}) - \text{mL}(\text{sample})] \times \text{N}(\text{KOH}) \times 56.1}{\text{Wt of sample in gms}} + (\text{Acid Value})$$

Size of Sample:

Unless specifically mentioned, the following size of sample is recommended:

Hydroxyl Value Range	Weight of Sample In Grams
300 - 350	0.75
250 - 300	1.0
200 -250	1.25
150 - 200	1.5
100 - 150	2.0
50 -100	3.0
20 -50	5.0
0 - 20	10.0

References:

Extracted from USP23/NF18 (401), 1995.